S1

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18/602reb_b ms002

eb_b_oslonki_piersi_ms002.indd 1 09.08.2023 14:25

GB • Instructions for use: Keep nipple protectors clean. Before the first and each subsequent use, wash them in warm water with a mild detergent and rinse thoroughly. Disinfect nipple protectors in boiling water for about 3 minutes. Do not wash nipple protectors in the dishwasher. Do not use antibacterial agents for cleaning. Avoid rough, surface-scratching tools during use and cleaning. Nipple protectors can be disinfected in the attached case according to the instructions. Do not leave nipple protectors in direct sunlight or near a source of heat. Apply the protector only on a dry and clean breast. Do not use any creams or Vaseline. Check the product before each use. Throw away at the first signs of damage or weakness. Do not enlarge the holes and do not create additional holes in nipple protectors, as it may cause damage. Before using the product, it is recommended to consult a midwife or paediatrician. Warning! This product is not a toy. Keep out of reach of children. Clean and dry nipple protectors store in the case. Keep the instructions of use as it contains important information. Batch code is on the packaging. **Instruction on how to disinfect nipple protectors in the disinfection case NOTE! Failure to follow this instruction may result in burns or damage to the product and/or the microwave oven.** • Always check the water level in the disinfection case before use.• Two nipple protectors in the disinfection case can be disinfected at one time. • This product is not a toy. Use and store the product out of the reach of children. **USE** 1. Before the first and each subsequent use, wash the disinfection case in warm water with a mild detergent, rinse thoroughly. 2. Pour 25 ml of boiled water into the base. ***IMPORTANT Adding too little water may damage the case. Adding too much water may result in inappropriate disinfection.*** 3. Put nipple protectors in the disinfection case. 4. Put the lid on and close tightly. 5. Place the closed case in the microwave oven. 6. For **750-1000 W** microwave ovens, the disinfection time is **2 minutes**. Do not exceed the disinfection time. **NOTE!** • The case must be placed on an even surface during disinfection, • Always wait at least 5 minutes for the disinfection case to cool down before removing it

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